

Product Information

Wet-End



INVADERM® VP – Dyeing auxiliary for dye penetration and levelling

Edition 2018

Features

INVADERM® VP is a dyeing auxiliary with affinity for the dye, which provides outstanding penetration and dye levelling effects.

Lightfastness: No detrimental effect on finished leathers.

Yellowing: No detrimental effect on finished leathers.

Stability to electrolytes: Very good.

Storage stability: Approx. 1 year if protected from frost and heat above 35°C.

Benefits

Owing to its retarding effect on dyes, **INVADERM® VP** improves the compatibility and penetration of the dyes, thus considerably improving both levelness and reproducibility.

With its affinity for the dye, **INVADERM® VP** improves the buffability of nubuck and suede as it ensures a more uniform penetration.

Application

In general, 0.5 – 2.0% **INVADERM® VP** (calculated on the shaved weight) is applied 10 minutes prior to the dye. However, for the very best effect on dye compatibility and penetration, **INVADERM® VP** can be added with the pre-dissolved dye.

Trade form at 20°C

turbid liquid

Colour at 20°C

yellowish to brown

Dry Content (%)

min. 53

pH – value (5%)

approx. 6

This information is for guidance only.
A Product Specification is available on request.

Regarding toxicological and ecological data see the safety data sheet of above product.

INVADERM® is a registered or filed trademark owned by or licensed to TFL in most countries.

Our application recommendations are in line with our present state of knowledge. They do not, however, exempt the customer from performing his own tests to determine the suitability of the supplied products for their intended purposes. Application of the products lies outside the scope of our control and therefore comes within the customer's sphere of responsibility. We guarantee the satisfactory quality of our products subject to general terms of sale and delivery.